# Galaxy Watch6 Classic



## Classic style. Latest innovation.

Meet the stylish smartwatch that's packed with the latest innovation — Samsung Galaxy Watch6 Classic. With a sleek silhouette, a rotating bezel and our largest display yet,<sup>1</sup> it boasts the iconic look of a traditional watch. Take your style on the go and use your Watch to call, text, pay and check notifications anytime, anywhere.<sup>2</sup> Keep moving and motivated with wellness insights right on your wrist, and develop better sleep habits with Advanced Sleep Coaching.<sup>3</sup>

The style you want. The innovation you need. That's Galaxy Watch6 Classic.

## The return of iconic style

Boost your style with a classic look — our latest smartwatch features a rotating bezel that you can turn to easily access your apps and more.

## Put your style on display

Show off your style with our largest display yet. Stunning and tough, this eyecatching display is made from durable crystal glass.

## A sleep coach on your wrist

Advanced Sleep Coaching<sup>3</sup> helps you sleep better so you wake up ready to take on your day. Use your Watch to plan your bedtime, detect snoring<sup>4</sup> and understand and track your sleep stages (awake, light, deep, REM).

#### Keep your goals on track

Work out to yourfull potential with fitness tracking<sup>3</sup> on Galaxy Watch6 Classic. Get useful performance insights on duration, distance, calories burned and more.

## Always-on heart monitoring

Get ahead of your heart health. Your Watch continually scans your heart rate to inform you when it's detected an irregular rhythm.

#### Stay connected on the go

Your Watch syncs with your other Galaxy devices seamlessly. So, have the power to call, text, take pics and play music all from your wrist.<sup>2</sup>

## SAMSUNG





## Galaxy Watch6 Classic

Available in LTE and Bluetooth





#### Galaxy Watch6 Classic

Available in LTE and Bluetooth

Colors	Case Black Silver Band Black Silver	Case Black Silver Band Black Silver
Dimensions	46.5mm x 46.5mm x 10.9mm <sup>6</sup>	42.5mm x 42.5mm x 10.9mm <sup>6</sup>
Weight	85g <sup>7</sup>	77g <sup>7</sup>
Display	1.47" (480px x 480px), 327ppi, Sapphire Crystal Glass	1.31" (432px x 432px), 330ppi, Sapphire Crystal Glass
Durability	5ATM / IP68 / MIL-STD-810H	5ATM / IP68 / MIL-STD-810H
Material	Stainless Steel	Stainless Steel
Health Features <sup>3</sup> (Powered with BioActive Sensor)	Advanced Sleep Coaching ECG / Irregular Heart Rhythm Notification <sup>5</sup> Personalized Heart Rate Zones BIA (Bioelectrical Impedance Analysis) Skin Temperature / Cycle Tracking <sup>8</sup> Blood Oxygen	Advanced Sleep Coaching ECG / Irregular Heart Rhythm Notification <sup>5</sup> Personalized Heart Rate Zones BIA (Bioelectrical Impedance Analysis) Skin Temperature / Cycle Tracking <sup>8</sup> Blood Oxygen
Battery	425mAh	300mAh
Connectivity	BT 5.2 / Wi-Fi 2.4GHz and 5GHz / GPS / NFC / LTE <sup>9</sup>	BT 5.2 / Wi-Fi 2.4GHz and 5GHz / GPS / NFC / LTE <sup>9</sup>
OS	Wear OS powered by Samsung	Wear OS powered by Samsung
Memory (RAM) and Storage	2GB RAM, 16GB internal storage	2GB RAM, 16GB internal storage
Sensors	Accel, Gyro, Barometer, Ambient Light, Compass	Accel, Gyro, Barometer, Ambient Light, Compass
Band	Sports Band M / L Size: 20mm (Hybrid Eco-Leather)	Sports Band S / M Size: 20mm (Hybrid Eco-Leather)
Physical User Interface	2 button + physical rotating bezel	2 button + physical rotating bezel
Miscellaneous	Mic, speaker, motor	Mic, speaker, motor

47mm Galaxy Watch6 Classic model only. <sup>2</sup>Compatibility for Galaxy wearable apps: Android 8.0 or later, 1.5GB or RAM. <sup>3</sup>Requires Samsung Health app version 6.22 or later. <sup>4</sup>Phone must be within 2 feet of user with some detection enabled through Samsung Health app version 6.22 or later. <sup>4</sup>Requires Samsung Galaxy Smartphone with Android 7.0 or later, and Samsung Health Monitor app (available only at the Samsung Galaxy app store). The ECG Monitor App intended for or adults 22 years and older analyzes pulse rate data to identify jepisodes of irregular heart rhythms suggestive of AFib and provides a notification is not intended to project of AFib and the absence of a notification is not intended to project on every episode of irregular rhythm suggestive of AFib and the absence of a notification is not intended to indicate no disease processis present; rather the feature is intended to opportunistically acquire pulse rate data when the user is still and analyze the data when determined sufficient toward surfacing anotification. The ECG Monitor App is not intended to replace traditional methods of diagnosis or treatments. The ECG app is not intended to replace traditional methods of diagnosis or treatments in the ECG app is not intended to replace traditional methods of diagnosis or treatments the ECG app is not intended for or users with known arrhythmias other than AFib or users under 22 years old. Users should not intended to replace 22 years old. Users should not intended to replace 22 years old. Users should not intended to replace 22 years old. Users should not intended to replace 22 years old. Users should not intended to replace 22 years old. Users should not intended to replace 22 years old. Users should not intended to replace 22 years old. Users should not replace 22 years o